
Compatibility Matrix – 2544 Cordless Multiport Expander

• Works with up to 5 wireless devices • Operting range 8 – 12 metres • 6.35mm N/O plug for conection to N/C system • 6.35mm socket for connecting a nurse call pendant • RJ10 socket for connecting hardwired pads & mats • 3 x AAA Battreies • 2406M 6V Plug Pack (optional) • Easy to program • Low Battery and Pad Lost indication • Available in Tip and Ring configuration – 2544TR • Available in Ring and Sleeve configuration – 2544RS		Self Reset	**2544 Cordless Multiport Expander** Suitable for nursing homes and hospitals where nurse call systems are installed. **What's in the box:** 1 x 2544 Multiport Expander 1 x Quick Start Guide 3 x High quality AAA batteries	
				12 Month Warranty
3515	Cordless Under Mattress Bed Pad 1 Year (830 x 260)	✓		Place **UNDER THE MATTRESS** Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 2544 Multiport Expander ****Requires 2719 Transmitter to function****
2629	Cordless Bed Pad 1 Year (760 x 500)	✓		Place **ON TOP OF THE MATTRESS**. Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 2544 Multiport Expander. ****Requires 2719 Transmitter to function****
2630	Cordless Chair Pad 1 Year (380 x 250)	✓		Sends **OUT OF CHAIR** message when getting out of chair. Resets **OUT OF CHAIR** when back on chair or via 2544 Multiport Expander. ****Requires 2719 Transmitter to function****
2624	Cordless Floor Mat 1 Year (1220 x 620)			Position in an area of **FALLS RISK**. Initiates alarm when **stepped on**. Reset the alarm via 2544 Multiport Expander. ****Requires 2719 Transmitter to function****
2719	Cordless Pad Transmitter			Required by Cordless Pads / Mats. Sends a signal to the 2544 Multiport Expander. ****1 year battery life. Replace the transmitter after one year****
2690	Universal Single Fold Bed Pad 18 Months (760 x 255)	✓		Place **ON TOP OF THE MATTRESS** Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 2544 Multiport Expander. ****Requires 2720 Transmitter to function****

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.

2691	Universal Chair Pad 18 Months (380 x 255)	✓		Sends **OUT OF CHAIR** message when getting out of chair. Resets **OUT OF CHAIR** when back on chair or via 2544 Multiport Expander. ****Requires 2720 Transmitter to function****
2689	Universal Under Mattress Bed Pad 18 Months (830 x 260)	✓		Place **UNDER THE MATTRESS** Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 2544 Multiport Expander. ****Requires 2720 Transmitter to function****
2692	Universal Floor Mat			Positioned strategically in an area of risk. Initiates alarm when stepped on. Reset at nurse call point. *Requires 2720 Transmitter to function*
2720	Universal Cordless Pad Transmitter			Required by Universal Pads/Mats if not connected via lead. Sends a signal to the 2544 Multiport Expander. *18 month battery life*
2670	Cordless PreAlert Bed Pad 1 Year (830 x 260)	✓		Located on **TOP OF THE MATTRESS**. Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 2544 Multiport Expander. ****Transmitter Included****
2671	Universal Single Fold Floor Mat 18 Months (900 x 600)			Position in an area of **FALLS RISK**. Initiates alarm when stepped on. Reset on 2544 Multiport Expander. Can be folded for storage or transportation. ****Requires 2720 Transmitter to function for cordless operation****
2672	Universal Dual Fold Floor Mat 18 Months (1190 x 600)			Position in an area of **FALLS RISK**. Initiates alarm when stepped on. Reset on 2544 Multiport Expander. Can be folded for storage or transportation. ****Requires 2720 Transmitter to function for cordless operation****
2628	Premium Chair Pad 1 Year (380 x 250)	✓		Position **ON CHAIR / WHEELCHAIR**. Sends **OUT OF CHAIR** message when getting out of chair. Resets **OUT OF CHAIR** when back on chair or via 2544 Multiport Expander.

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.

2626	Premium Bed Pad 1 Year (760 x 500)	✓		Position **ABOVE MATTRESS** Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 2544 Multiport Expander.
3516	Premium Under Mattress Bed Pad 1 Year (830 x 260)	✓		Position **UNDER MATTRESS** Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 2544 Multiport Expander.
2633	Premium Heavy-Duty Floor Mat 1 Year (1220 x 610)			Position in an area of **falls risk**. The 2544 Multiport Expander alarms when stood on. Reset via 2544 Multiport Expander.
2588	Personal Emergency Transmitter (PET)			Initiates alarm when **RED Button** is pressed. Reset on 2544 Multiport Expander.
2721	iBeam			Sends alarm when triggered. Features a narrow detection zone to reduce false alarms, triggers when a patient breaks the beam. Reset on 2544 Multiport Expander.
2638	Nurse Call Pendant			Simple operation. Plugs into the **2544** and **2549** Multiport expanders. **N / O** push button. 2 metre flexible cord.
2406M	cura1 6V Plug Pack			Plug Pack with 1.3mm connection tip.
2558	Wireless Waterproof Call Button w/ pullcord	✓		Initiates alarm when **RED button** is pressed or **cord is pulled**. Includes wireless **reset button** or reset on **2544 Multiport Expander**
2595	Wristband Alarm Transmitter			Staff or patients wear on the wrist. Use as a duress alarm or call button. Initiates alarm when **RED Button** is pressed. Reset on 2544 Multiport Expander.

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.

3502	ActiveCare Wireless Handset	✓		Initiates alarm when **RED button** is pressed. Includes **GREEN** wireless reset button or reset via 2544 Multiport Expander.
3503	Wireless Door Contact			Sends alarm when door is **opened**. Can also be used as a simple **push button** call point. Reset on 2544 Multiport Expander.
2698	Universal Transmitter	✓		Turn any third-party alarm device into a wireless transmitter. Can also be used as a wireless **door contact**. Togglable **Reset** button.
2599	Wireless SOS Call Button			Wall mountable emergency call point. Initiates alarm when **RED Button** is pressed. Reset via 2544 Multiport Expander.
2730	SnapOn Bed Pad 1 Year (800 x 255)			Place **ON TOP OF THE MATTRESS** Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via connected monitor. ****Requires 2732 SnapOn Pad Transmitter or 2733 SnapOn Pad Cable to function****
2731	SnapOn Chair Pad 1 Year (420 x 180)			Sends **OUT OF CHAIR** message when getting out of bed. Resets **OUT OF CHAIR** when back in bed or via connected monitor. ****Requires 2732 SnapOn Pad Transmitter or 2733 SnapOn Pad Cable to function****
2732	SnapOn Pad Transmitter			Required by SnapOn Pads. Sends a signal to the 2544 Multiport Expander. ****Replaceable CR2045 battery.****
2733	SnapOn Pad Cable 2m			Easily connects to the 2730 and 2731 SnapOn Pad range. Sends a signal to wired monitors through an RJ-10 connetor.
2669	Cordless Under Mattress PreAlert Bed Pad 1 Year (820 x 240)			Place **UNDER THE MATTRESS** Sends **OUT OF BED** message when at edge of bed. Resets **OUT OF BED** when back in bed or via 2700 Falls Monitor

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.