
• Durable ABS plastic enclosure with soft touch removable protection cover and Caregiver Reset Key • Amber LED visual alarm indicator • Caregiver secure reset option • Battery operated or optional external plug pack • 6 – 12 months battery life (nominal) • Easy to access battery compartment for staff • Easily connected to a nurse call system, selectable alarm tones, 0, 2 sec, and 15min selectable alarm delay settings • Easy thumbwheel adjustable alarm volume • "Pause" feature allows caregiver to attend to the patient without initiating an alarm (will auto reset after 3 minutes). • Alerts for lost nurse call connection, lost AC power & low battery • Change Pad alarm feature when SafeTime pads used		Self Reset	**2701 Falls Monitor with SafeTime** Suitable for private residences, hospitals and care homes, its versatility turns a basic room into a highly optioned care and monitoring facility offering multiple solutions in one simple to operate unit. **What's in the box:** 1 x 2701 Standard Falls Monitor 1 x Quick Start Guide 1 x 2642 Caregiver Reset Key 1 x Durable Soft Protective Cover 3 x Alkaline AA Batteries	
				24 Month Warranty
3516	Premium Under Mattress Bed Pad 1 Year (830 x 260)	✓		Position **UNDER MATTRESS** Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 2701 Falls Monitor.
2626	Premium Bed Pad 1 Year (760 x 500)	✓		Position **ABOVE MATTRESS** Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 2701 Falls Monitor.
2628	Premium Chair Pad 1 Year (380 x 250)	✓		Position **ON CHAIR / WHEELCHAIR**. Sends **OUT OF CHAIR** message when getting out of chair. Resets **OUT OF CHAIR** when back on chair or via 2701 Falls Monitor.
2633	Premium Heavy-Duty Floor Mat 1 Year (1220 x 610)			Position in an area of **falls risk**. The 2701 Falls Monitor alarms when stood on. Reset via 2701 Falls Monitor.
2622	SafeTime Bed Pad 90 Days (760 x 500)	✓		Position **ABOVE MATTRESS** Sends **OUT OF BED** message when getting out of bed. Alarms when the pad **expires**. Resets **OUT OF BED** when back in bed or via 2701 Falls Monitor.
2620	SafeTime Wheelchair / Chair Pad - 90 days (280 x 250)	✓		Position **ON CHAIR / WHEELCHAIR**. Sends **OUT OF CHAIR** message when getting out of chair. Alarms when the pad **expires**. Resets **OUT OF CHAIR** when back on chair or via 2701 Falls Monitor.
2690	Universal Single Fold Bed Pad 18 Months (760 x 255)	✓		Place **ON TOP OF THE MATTRESS** Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 2701 Falls Monitor. ****Requires 2693 USB Break Cord to function****

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.

2691	Universal Chair Pad 18 Months (380 x 255)	✓		Sends **OUT OF CHAIR** message when getting out of chair. Resets **OUT OF CHAIR** when back on chair or via 2701 Falls Monitor. ****Requires 2693 USB Break Cord to function****
2689	Universal Under Mattress Bed Pad 18 Months (830 x 260)	✓		Place **UNDER THE MATTRESS** Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 2701 Falls Monitor. ****Requires 2693 USB Break Cord to function****
2671	Universal Single Fold Floor Mat 18 Months (900 x 600)			Position in an area of **FALLS RISK**. Initiates alarm when stepped on. Reset on 2701 Falls Monitor. Can be folded for storage or transportation. ****Requires 2693 USB Break Cord to function****
2672	Universal Dual Fold Floor Mat 18 Months (1190 x 600)			Position in an area of **FALLS RISK**. Initiates alarm when stepped on. Reset on 2701 Falls Monitor. Can be folded for storage or transportation. ****Requires 2693 USB Break Cord to function****
2692	Universal Floor Mat 18 Months (1220 x 610)			Position in an area of **FALLS RISK**. Initiates alarm when stepped on. Reset the alarm via 2701 Falls Monitor. ****Requires 2693 USB Break Cord to function****
2742	Spare Soft Cover for Monitors			Soft touch removable protection case provides extra protection for the 2701 Monitor.
2406	cura1 6V Plug Pack			Optional **6 Volt DC** plug pack.
2642	Caregiver Reset Key			Allows caregivers to **reset alarms**. Allows authorised staff to **program** various cura1 equipment including **3500R Monitors** Key chain can be **attached** to belt or clothing. **Lightweight**.
3504	Mounting Bracket			2701 Falls Monitor slides into bracket. Bracket mounts onto surfaces with **double-sided tape** or **M3 screws**.

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.

2636	Easy Release Wheelchair / Chair Seatbelt	✓		Activates an alarm if buckle is released.
2643	Nurse Call Cable 6.35 Mono to 6.35 Mono 2m			Easily connected to a nurse call system. (Some nurse call systems may require a different cable. Please contact customer support for assistance.)
2646	Nurse Call Cable 6.35 Mono to 6.35 Mono 3m			Easily connected to a nurse call system. (Some nurse call systems may require a different cable. Please contact customer support for assistance.)
2493	Nurse Call Cable 6.35 Mono to 3.5 Mono 2m			Easily connected to a nurse call system. (Some nurse call systems may require a different cable. Please contact customer support for assistance.)
2461	Nurse Call Cable 6.35 Mono to 6.35 T&S 2m			Easily connected to a nurse call system. (Some nurse call systems may require a different cable. Please contact customer support for assistance.)
2460	Nurse Call Cable 6.35 T&S to 6.35 T&S 2m			Easily connected to a nurse call system. (Some nurse call systems may require a different cable. Please contact customer support for assistance.)
2730	SnapOn Bed Pad 1 Year (800 x 255)			Place **ON TOP OF THE MATTRESS** Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via connected monitor. ****Requires 2733 SnapOn Pad Cable to function****
2731	SnapOn Chair Pad 1 Year (420 x 180)			Sends **OUT OF CHAIR** message when getting out of bed. Resets **OUT OF CHAIR** when back in bed or via connected monitor. ****Requires 2733 SnapOn Pad Cable to function****
2733	SnapOn Pad Cable 2m			Easily connects to the 2730 and 2731 SnapOn Pad range. Sends a signal to wired monitors through an RJ-10 connetor.

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.