
	• Works with up to 5 wireless devices • Operating range 8 – 12 metres • RJ12 socket for connecting to N/C system • 3 x AAA Batteries • 2406M 6V Plug Pack (optional) • Easy to program • Low Battery and Pad Lost indication	Self Reset	**2549 Cordless Multiport Expander with Bracket Mount** Suitable for nursing homes and private/public hospitals where nurse call systems are installed. **What's in the box:** 1 x 2549 Multiport Expander Bracket Mount 1 x 3504 Mounting Bracket 1 x Quick Start Guide 3 x AAA Batteries 1 x Lead to connect to N/C system	
				12 Month Warranty
3515	Cordless Under Mattress Bed Pad 1 Year (830 x 260)	✓		Place **UNDER THE MATTRESS** Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 2549 Multiport Expander ****Requires 2719 Transmitter to function****
2629	Cordless Bed Pad 1 Year (760 x 500)	✓		Place **ON TOP OF THE MATTRESS**. Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 2549 Multiport Expander. ****Requires 2719 Transmitter to function****
2630	Cordless Chair Pad 1 Year (380 x 250)	✓		Sends **OUT OF CHAIR** message when getting out of chair. Resets **OUT OF CHAIR** when back on chair or via 2549 Multiport Expander. ****Requires 2719 Transmitter to function****
2624	Cordless Floor Mat 1 Year (1220 x 620)			Position in an area of **FALLS RISK**. Initiates alarm when **stepped on**. Reset the alarm via 2549 Multiport Expander. ****Requires 2719 Transmitter to function****
2719	Cordless Pad Transmitter			Required by Cordless Pads / Mats. Sends a signal to the 2549 Multiport Expander. ****1 year battery life. Replace the transmitter after one year****

CM016-7

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.

2690	Universal Single Fold Bed Pad 18 Months (760 x 255)	✓		Place **ON TOP OF THE MATTRESS** Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 2549 Multiport Expander. ****Requires 2720 Transmitter to function****
2689	Universal Under Mattress Bed Pad 18 Months (830 x 260)	✓		Place **UNDER THE MATTRESS** Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 2549 Multiport Expander. ****Requires 2720 Transmitter to function****
2691	Universal Chair Pad 18 Months (380 x 255)	✓		Sends **OUT OF CHAIR** message when getting out of chair. Resets **OUT OF CHAIR** when back on chair or via 2549 Multiport Expander. ****Requires 2720 Transmitter to function****
2692	Universal Floor Mat 18 Months (1220 x 610)			Position in an area of **FALLS RISK**. Initiates alarm when stepped on. Reset the alarm via 2549 Multiport Expanderr. ****Requires 2720 Transmitter to function****
2720	Universal Cordless Pad Transmitter			Required by Universal Pads / Mats if not connected via lead. Sends a signal to the 2549 Multiport Expander. ****18 month battery life. Replace the transmitter after 18 months****
2670	Cordless PreAlert Bed Pad 1 Year (830 x 260)	✓		Located on **TOP OF THE MATTRESS**. Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 2549 Multiport Expander. ****Transmitter Included****
2671	Universal Single Fold Floor Mat 18 Months (900 x 600)			Position in an area of **FALLS RISK**. Initiates alarm when stepped on. Reset on 2549 Multiport Expander. Can be folded for storage or transportation. ****Requires 2720 Transmitter to function for cordless operation****

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.

2558	Wireless Waterproof Call Button w/pullcord	✓		Initiates alarm when **RED button** is pressed or **cord is pulled**. Includes wireless **reset button** or reset on **2549 Multiport Expander**
2595	Wristband Alarm Transmitter			Staff or patients wear on the wrist. Use as a duress alarm or call button. Initiates alarm when **RED Button** is pressed. Reset on 2549 Multiport Expander.
2599	Wireless SOS Call Button			Wall mountable emergency call point. Initiates alarm when **RED Button** is pressed. Reset via 2549 Multiport Expander.
3502	ActiveCare Wireless Handset	✓		Initiates alarm when **RED button** is pressed. Includes **GREEN** wireless reset button or reset via 2549 Multiport Expander.
3503	Wireless Door Contact			Sends alarm when door is **opened**. Can also be used as a simple **push button** call point. Reset on 2549 Multiport Expander.
2698	Universal Transmitter	✓		Turn any third-party alarm device into a wireless transmitter. Can also be used as a wireless **door contact**. Togglable **Reset** button.
3504	Mounting Bracket			2549 Multiport Expander slides into bracket. Bracket mounts onto surfaces with **double-sided tape** or **M3 screws**.
2669	Cordless Under Mattress PreAlert Bed Pad 1 Year (820 x 240)			Place **UNDER THE MATTRESS** Sends **OUT OF BED** message when at edge of bed. Resets **OUT OF BED** when back in

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.

2672	Universal Dual Fold Floor Mat 18 Months (1190 x 600)			Position in an area of **FALLS RISK**. Initiates alarm when stepped on. Reset on 2549 Multiport Expander. Can be folded for storage or transportation. ****Requires 2720 Transmitter to function for cordless operation****
2588	Personal Emergency Transmitter (PET)			Initiates alarm when **RED Button** is pressed. Reset on 2549 Multiport Expander.
2721	iBeam			Sends alarm when triggered. Features a narrow detection zone to reduce false alarms, triggers when a patient breaks the beam. Reset on 2549 Multiport Expander.
2638	Nurse Call Pendant			Simple operation. Plugs into the **2544** and **2549** Multiport expanders. **N/O** push button. 2 metre flexible cord.
2406M	cura1 6V Plug Pack			Plug Pack with 1.3mm connection tip.
4160	Curly Cord RJ12 to 6.35 R/A Plug St 300mm T&S			Connects the 2549 directly to a nurse call system.
4161	Curly Cord RJ12 to 6.35 R/A Plug St 300mm T&R			Connects the 2549 directly to a nurse call system.
4162	Curly Cord RJ12 to 6.35 R/A Plug St 300mm R&S			Connects the 2549 directly to a nurse call system.
4163	Curly Cord RJ12 to 3.5 R/A Plug St 300mm T&S			Connects the 2549 directly to a nurse call system.
4164	Curly Cord RJ12 to 3.5 R/A Plug St 300mm T&R			Connects the 2549 directly to a nurse call system.
4165	Curly Cord RJ12 to 3.5 R/A Plug St 300mm R&S			Connects the 2549 directly to a nurse call system.

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.

2730	SnapOn Bed Pad 1 Year (800 x 255)			Place **ON TOP OF THE MATTRESS** Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via connected monitor. ****Requires 2732 SnapOn Pad Transmitter to function****
2731	SnapOn Chair Pad 1 Year (420 x 180)			Sends **OUT OF CHAIR** message when getting out of bed. Resets **OUT OF CHAIR** when back in bed or via connected monitor. ****Requires 2732 SnapOn Pad Transmitter to function****
2732	SnapOn Pad Transmitter			Required by SnapOn Pads. Sends a signal to the Multiport Expander. ****Replaceable CR2045 battery.****

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.