
		Self Reset	**3010 SeizureAlert Monitor**	
• Detect & raise an alarm as a seizure occurs • 6.35mm stereo T & S output for connecting to a nurse call system • Battery backup built in • Easy to program			Suitable for nursing homes and private/public hospitals where nurse call systems are installed or a home environment. **What's in the box:** 1 x 3010 SeizureAlert Monitor 1 x Instruction manual 1 x Nurse Call cable 1 x Plug pack	
				12 Month Warranty
3508	Activecare LED Pager			Reports alerts from up to **3** wireless transmitters. Can be mounted on a wall or carried by a caregiver. Alarms via **sound** or **vibrate** and **Red LED** illuminates.
3515	Cordless Under Mattress Bed Pad 1 Year (830 x 260)	✓		Place **UNDER THE MATTRESS** Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 3508 LED Pager ****Requires 2719 Transmitter to function****
2629	Cordless Bed Pad 1 Year (760 x 500)	✓		Place **ON TOP OF THE MATTRESS**. Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 3508 LED Pager. ****Requires 2719 Transmitter to function****
2630	Cordless Chair Pad 1 Year (380 x 250)	✓		Sends **OUT OF CHAIR** message when getting out of chair. Resets **OUT OF CHAIR** when back on chair or via 3508 LED Pager. ****Requires 2719 Transmitter to function****
2624	Cordless Floor Mat 1 Year (1220 x 620)			Position in an area of **FALLS RISK**. Initiates alarm when **stepped on**. Reset the alarm via 3508 LED Pager. ****Requires 2719 Transmitter to function****
2588	PET (Personal Emergency Transmitter)			Initiates alarm when **RED Button** is pressed. Reset on 3508 LED Pager.
2719	Cordless Pad Transmitter			Required by Cordless Pads / Mats. Sends a signal to the 3508 LED Pager. ****1 year battery life. Replace the transmitter after one year****

CM013-7

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.

2690	Universal Single Fold Bed Pad 18 Months (760 x 255)	✓		Place **ON TOP OF THE MATTRESS** Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 2700 Falls Monitor. ****Requires 2720 Transmitter to function****
2689	Universal Under Mattress Bed Pad 18 Months (830 x 260)	✓		Place **UNDER THE MATTRESS** Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 2700 Falls Monitor. ****Requires 2720 Transmitter to function****
2691	Universal Chair Pad 18 Months (380 x 255)	✓		Sends **OUT OF CHAIR** message when getting out of chair. Resets **OUT OF CHAIR** when back on chair or via 2700 Falls Monitor. ****Requires 2720 Transmitter to function****
2692	Universal Floor Mat 18 Months (1220 x 610)			Position in an area of **FALLS RISK**. Initiates alarm when stepped on. Reset the alarm via 2700 Falls Monitor. ****Requires 2720 Transmitter to function****
2720	Universal Cordless Pad Transmitter			Required by Universal Pads / Mats if not connected via lead. Sends a signal to the 2700 Falls Monitor. ****18 month battery life. Replace the transmitter after 18 months****
2670	Cordless PreAlert Bed Pad 1 Year (830 x 260)	✓		Located on **TOP OF THE MATTRESS**. Sends **OUT OF BED** message when getting out of bed. Resets **OUT OF BED** when back in bed or via 2700 Falls Monitor. ****Transmitter Included****
2671	Universal Single Fold Floor Mat 18 Months (900 x 600)			Position in an area of **FALLS RISK**. Initiates alarm when stepped on. Reset on 2700 Falls Monitor. Can be folded for storage or transportation. ****Requires 2720 Transmitter to function for cordless operation****
2672	Universal Dual Fold Floor Mat 18 Months (1190 x 600)			Position in an area of **FALLS RISK**. Initiates alarm when stepped on. Reset on 2700 Falls Monitor. Can be folded for storage or transportation. ****Requires 2720 Transmitter to function for cordless operation****

All radio frequency (RF) signal transmission and reception distances quoted are indicative only and based on testing carried out under ideal conditions. Actual transmission and reception distances may vary and be dependent on specific environmental conditions.